

PRODUCT INFORMATION

Common Name	AMG 477, REMD-477
Synonyms	Glucagon receptor;GL-R
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	IgG2
Reactivity	Human
Target	GCGR
Uniprot ID	P47871
Description	Anti-GCGR(volagidemab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only. All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.
Conjugate	Unconjugated






Figure 1. Flow cytometry analysis with 1 μ g/mL Anti-GCGR (volagidemab biosimilar) mAb (BME100142) on Expi293 cells transfected with Human GCGR protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

