

PRODUCT INFORMATION

Common Name	hLIV22
Conjugate	Unconjugated
Synonyms	LIV1;SLC39A6;ZIP-6
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-10000
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	LIV-1
Uniprot ID	Q13433
Description	Anti-LIV-1(ladiratuzumab biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.



Anti-LIV-1 (ladiratuzumab biosimilar) mAb ELISA

0.2 μ g of Human LIV-1, His tagged protein per well




Figure 1. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human LIV-1 Protein, His Tag ([getskuurl sku="PME100759"]) can bind Anti-LIV-1 Neutralizing antibody (BME100113) in a linear range of 0.92-58.59 ng/mL.

DIMABIO CONFIDENTIAL

