

PRODUCT INFORMATION

Clone ID	9E10
Target	Myc
Synonyms	9E10 antibody; 9E10 anti myc; 9E10 myc; anti-Myc; anti-myc antibodies; Anti-myc Antibody; anti c-myc; anti cmc; anti myc; c-myc antibody; cmc antibody; myc antibody; tag antibody
Host Species	Mouse
Description	Anti-Myc Monoclonal Antibody
Delivery	In Stock
Uniprot ID	N/A
IgG type	Mouse IgG
Clonality	Monoclonal
Reactivity	Myc tag
Applications	FC, IF, IP, WB
Recommended Dilutions	WB: 1:500-1:2000, FC: 1:100, IF 1:100, IP 1:100-200
Purification	N/A
Formulation & Reconstitution	0.64 mg/ml in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
Storage & Shipping	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

