

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	O00206
<b>Common Name</b>	EB 05;EB05;EB-05;EB 5;EB-5;EB5;EB-07;EB07;EB-07;EB-7;EB7;EB 7;NI-0101;NI 0101;NI0101;NI-101;NI101;NI 101
<b>Synonyms</b>	TOLL;CD284;TLR-4;ARMD10
<b>Applications</b>	N/A
<b>Recommended Dilutions</b>	N/A
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Whole mAb
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	TLR4
<b>Description</b>	Anti-TLR4(paridiprubart biosimilar) IgG1 mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated

DIMA BIOTECH CONFIDENTIAL

