

**PRODUCT INFORMATION**

<b>Target</b>	CD46, CD73, CDH17, CEACAM5, EPCAM, FZD10, GPA33, GUCY2C, HER3
<b>Description</b>	The embedded and sliced tumor sections (4µm section) from Balb/C nu mice implanted with HT55 cells
<b>Host Cells</b>	HT55
<b>Applications</b>	IHC
<b>Host mice</b>	Balb/C nu
<b>Delivery</b>	1 week
<b>Guarantee</b>	36 months from the date of manufacture
<b>Background</b>	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
<b>Storage &amp; Shipping</b>	Room temperature
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

