

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P05067 |
| Common Name | LY-3372993;LY 3372993;LY3372993;N3pG 4;N3pG4;N3PG-IV |
| Synonyms | AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | APP |
| Description | Biotinylated Anti-APP(remternetug biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | Biotinylated |

DIMABIO CONFIDENTIAL

