

PRODUCT INFORMATION

Uniprot ID	P05113
Common Name	CDP-835;CDP835;CDP 835;CEP-38072;CEP38072;CEP 38072;CINQAERO;CINQAIR;Cinquil;CTx-55700;CTx55700;Ctx 55700;DCP-835;DCP835;DCP 835;SCH 55700;SCH55600;SCH-55700;TRFK-5;TRFK 5;TRFK5
Synonyms	EDF;TRF;IL-5;Interleukin-5
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	IL5
Description	Biotinylated Anti-IL5(reslizumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

