

PRODUCT INFORMATION

Uniprot ID	P08887
Common Name	Keyzara;REGN-88;REGN88;REGN 88;SAR-153191;SAR153191;SAR153191
Synonyms	IL6Q;gp80;CD126;HIES5;IL-6R;IL6RA;IL6RQ;IL-1Ra;IL-6RA;IL6QTL;IL-6R-1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	IL6R
Description	Biotinylated Anti-IL6R(sarilumab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Biotinylated

DIMABIO CONFIDENTIAL

