

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P10636 |
| Common Name | JNJ-63733657;JNJ 63733657;JNJ63733657;JNJ-3657;JNJ3657;JNJ 3657 |
| Synonyms | TAU;FTD1;MSTD;PPND;DDPAC;MAPTL;MTBT1;MTBT2;tau-40;FTDP-17;PPP1R103;Tau-PHF6 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | MAPT |
| Description | Biotinylated Anti-MAPT(posdinemab biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | Biotinylated |

DIMABIO CONFIDENTIAL

