

PRODUCT INFORMATION

Uniprot ID	Q9NZQ7
Common Name	STI-1014;STI1014;STI 1014;STI-A1014;STI A1014;ZKAB-001;ZKAB001;ZKAB 001;ZKAB-1;ZKAB1;ZKAB 1
Synonyms	B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	PD-L1
Description	Biotinylated Anti-PD-L1(socazolimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Biotinylated

