

**PRODUCT INFORMATION**

|   |   |
|---|---|
| <b>Clone ID</b>                         | 34B1  |
| <b>Target</b>                           | Rabbit IgG  |
| <b>Synonyms</b>                         | N/A   |
| <b>Host Species</b>                     | Rabbit  |
| <b>Description</b>                      | Biotinylated Anti-Rabbit IgG antibody(34B1), Goat mAb   |
| <b>Delivery</b>                         | 2-3 weeks   |
| <b>Uniprot ID</b>                       | N/A   |
| <b>IgG type</b>                         | Goat IgG  |
| <b>Clonality</b>                        | Monoclonal  |
| <b>Reactivity</b>                       | Rabbit  |
| <b>Applications</b>                     | ELISA   |
| <b>Recommended Dilutions</b>            | ELISA 1/5000-10000  |
| <b>Purification</b>                     | Purified from cell culture supernatant by affinity chromatography   |
| <b>Formulation &amp; Reconstitution</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.   |
| <b>Storage &amp; Shipping</b>           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| <b>Background</b>                       | N/A   |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | Biotinylated  |

