

PRODUCT INFORMATION

Clone ID	4E12
Target	SIGLEC15
Synonyms	CD33 antigen-like 3;SIGLEC-15;CD33L3;sialic acid-binding Ig-like lectin 15;Siglec15;Siglec-15
Host Species	Rabbit
Description	Biotinylated Anti-SIGLEC15 antibody(4E12), IgG1 Chimeric mAb
Delivery	2-3 weeks
Uniprot ID	Q6ZMC9
IgG type	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1/100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	SIGLEC15 (Sialic Acid Binding Ig Like Lectin 15) is a Protein Coding gene. Diseases associated with SIGLEC15 include Osteoporosis;Juvenile and Osteoporosis. Among its related pathways are Innate Immune System and RET signaling. An important paralog of this gene is SIGLEC1.
Usage	Research use only
Conjugate	Biotinylated

