

PRODUCT INFORMATION

Target	ADAM9
Synonyms	CORD9;MCMP;MDC9;Mltng
Description	Recombinant Cynomolgus ADAM9 protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	A0A2K5X4X8
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	ADAM9(Ala29-Gly698) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 74.9 kDa after removal of the signal peptide. The apparent molecular mass of cADAM9-His is approximately 55-100 kDa due to glycosylation.
Package	10 µg \$258.00 ; 50 µg \$969.00 ; 100 µg \$1421.00
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. The protein encoded by this gene interacts with SH3 domain-containing proteins, binds mitotic arrest deficient 2 beta protein, and is also involved in TPA-induced ectodomain shedding of membrane-anchored heparin-binding EGF-like growth factor. Several alternatively spliced transcript variants have been identified for this gene. [provided by RefSeq, Jul 2010]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Cynomolgus ADAM9 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

