


PRODUCT INFORMATION

Target	CD123
Description	Monoclonal Cell Line Derived from K562 Cells, Engineered for Stable Expression of Human CD123 Using Lentiviral Technology
Host Cells	K562
Uniprot ID	P26951
Applications	FACS Data
Growth media	RPMI-1640+10% FBS+1% P.S+1% Gln+2 ug/mL Puromycin
Package	5E6 Cells/mL
Suggested Control	SKU: BME100003
Warranty and Disclaimer	1. Please inspect cells upon receipt and report any issues promptly. 2. We offer one-time replacements for issues reported within a week of receipt. 3. User-induced issues are not eligible for free replacements. 4. We do not accept liability for damages resulting from cell use, storage, or loss. 5. Feedback received more than one month after receipt will not be processed.
Storage & Shipping	Cells are shipped using dry ice and require liquid nitrogen storage for long term preservation.
Synonyms	IL3R; IL3RA; IL-3Ra; IL-3R-alpha; IL3RAY; IL3RX; IL3RY; CD123 antigen; CD123; hIL3Ra; hIL-3Ra; MGC34174; IL-3 R alpha
Background	Interleukin 3 receptor alpha (low affinity) (IL3RA); also known as CD123 (Cluster of Differentiation 123) is a 70-kD glycoprotein member of the hematopoietin receptor superfamily. This protein associates with a beta subunit common to the receptors for IL-5 and granulocyte-macrophage colony-stimulating factor (GM-CSF) to form a high-affinity receptor for IL-3. The interleukin-3 receptor a chain (CD123) has been identified as a potential immunotherapeutic target because it is overexpressed in AML compared with normal hematopoietic stem cells.
Usage	For research use only.



Hu_CD123 K562 Cell Line



-  Human IgG
-  Anti-CD123 (talacotuzumab biosimilar) mAb (SKU: BME100003)

Figure 1. Flow cytometry analysis of human CD123 overexpression using Hu_CD123 K562 Cell Line (Cat. No. CEL100011) and Anti-CD123 (talacotuzumab biosimilar) mAb (Cat. No. BME100003)

