

PRODUCT INFORMATION

Target	AREG
Synonyms	AR;CRDGF;AREGB;SDGF
Description	Recombinant human AREG protein with C-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P15514
Expression Host	HEK293
Tag	C-Human Fc Tag
Molecular Characterization	AREG(Ser20-Lys187) hFc(Glu99-Ala330)
Molecular Weight	The protein has a predicted molecular mass of 44.9 kDa after removal of the signal peptide. The apparent molecular mass of AREG-hFc is approximately 53-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	The protein encoded by this gene is a member of the epidermal growth factor family. It is an autocrine growth factor as well as a mitogen for astrocytes, Schwann cells and fibroblasts. It is related to epidermal growth factor (EGF) and transforming growth factor alpha (TGF-alpha). The protein interacts with the EGF/TGF-alpha receptor to promote the growth of normal epithelial cells, and it inhibits the growth of certain aggressive carcinoma cell lines. It also functions in mammary gland, oocyte and bone tissue development. This gene is associated with a psoriasis-like skin phenotype, and is also associated with other pathological disorders, including various types of cancers and inflammatory conditions. [provided by RefSeq, Apr 2014]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human AREG Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

