

PRODUCT INFORMATION

Target	CCL24
Synonyms	C-C motif chemokine 24 (CCL24);CK-beta-6;Eosinophil chemotactic protein;Eotaxin-2;Myeloid progenitor inhibitory factor 2;Small-inducible cytokine A24;Ckb-6;MPIF-2;SCYA24
Description	Recombinant Human C-C Motif Chemokine 24 is produced by our Mammalian expression system and the target gene encoding Val27-Cys119 is expressed with a 6His tag at the C-terminus.
Delivery	In Stock
Uniprot ID	O00175
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	Not available
Molecular Weight	11.5 KDa
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Formulation & Reconstitution	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Cytokine are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteine. CCL24 displays the chemotactic for resting T-lymphocytes and eosinophils. CCL24 has lower chemotatic activity for neutrophils, none for monocytes and activated lymphocytes. CCL24 is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell lines.
Usage	Research use only






Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

DIMABIO CONFIDENTIAL

