

PRODUCT INFORMATION

Target	CD72
Synonyms	Lyb-2
Description	Recombinant human CD72 protein with N-terminal human Fc tag
Delivery	In Stock
Uniprot ID	P21854
Expression Host	HEK293
Tag	N-Human Fc Tag
Molecular Characterization	hFc(Glu99-Ala330) CD72(Arg117-Asp359)
Molecular Weight	The protein has a predicted molecular mass of 54.2 kDa after removal of the signal peptide. The apparent molecular mass of hFc-CD72 is approximately 55-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Plays a role in B-cell proliferation and differentiation.[UniProtKB/Swiss-Prot Function]
Usage	Research use only






Figure 1. Human CD72 Protein, hFc Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

