

PRODUCT INFORMATION

Target	ELAPOR1
Synonyms	ELAPOR1;ndosome-Lysosome Associated Apoptosis And Autophagy Regulator 1;IG121;IAA1324
Description	Recombinant Human Endosome/Lysosome-associated Apoptosis and Autophagy Regulator 1 is produced by our Mammalian expression system and the target gene encoding Thr42-Lys910 is expressed with a 6His tag at the C-terminus.
Delivery	In Stock
Uniprot ID	Q6UXG2
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	Not available
Molecular Weight	96.2 KDa
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Formulation & Reconstitution	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Endosome/lysosome-associated apoptosis and autophagy regulator (ELAPOR1), also known as EIG121 protein, is a type I transmembrane protein induced by estrogen. The estrogen-induced gene 121 (EIG121) has been associated with breast and endometrial cancers, but its mechanism of action remains unknown. May protect cells from cell death by inducing cytosolic vacuolization and upregulating the autophagy pathway. That EIG121 is a good endometrial biomarker associated with a hyperestrogenic state and estrogen-related type I endometrial adenocarcinoma.
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Greater than 95% as determined by reducing SDS-PAGE.

DIMABIO CONFIDENTIAL

