

**PRODUCT INFORMATION**

<b>Target</b>	MERTK
<b>Synonyms</b>	MER;RP38;c-Eyk;c-mer;Tyro12
<b>Description</b>	Recombinant human MERTK Protein with C-terminal 10×His tag
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	Q12866
<b>Expression Host</b>	HEK293
<b>Tag</b>	C-10×His Tag
<b>Molecular Characterization</b>	MERTK(Ala21-Ile505) 10×His tag
<b>Molecular Weight</b>	The protein has a predicted molecular mass of 54.0 kDa after removal of the signal peptide.
<b>Purity</b>	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Background</b>	This gene is a member of the MER/AXL/TYRO3 receptor kinase family and encodes a transmembrane protein with two fibronectin type-III domains, two Ig-like C2-type (immunoglobulin-like) domains, and one tyrosine kinase domain. Mutations in this gene have been associated with disruption of the retinal pigment epithelium (RPE) phagocytosis pathway and onset of autosomal recessive retinitis pigmentosa (RP). [provided by RefSeq, Jul 2008]
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated






Figure 1. Human MERTK Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

