

**PRODUCT INFORMATION**

|   |   |
|---|---|
| <b>Target</b>                           | PLXDC2  |
| <b>Synonyms</b>                         | TEM7R   |
| <b>Description</b>                      | Recombinant Human PLXDC2 Protein with C-terminal 6×His tag  |
| <b>Delivery</b>                         | In Stock  |
| <b>Uniprot ID</b>                       | Q6UX71  |
| <b>Expression Host</b>                  | HEK293  |
| <b>Tag</b>                              | C-6×His Tag   |
| <b>Molecular Characterization</b>       | PLXDC2(Lys31-His453) 6×His tag  |
| <b>Molecular Weight</b>                 | The protein has a predicted molecular mass of 48.7 kDa after removal of the signal peptide. The apparent molecular mass of PLXDC2-His is approximately 55-100 kDa due to glycosylation.   |
| <b>Purity</b>                           | The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.  |
| <b>Formulation &amp; Reconstitution</b> | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.   |
| <b>Storage &amp; Shipping</b>           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| <b>Background</b>                       | May play a role in tumor angiogenesis.[UniProtKB/Swiss-Prot Function]   |
| <b>Usage</b>                            | Research use only   |
| <b>Conjugate</b>                        | Unconjugated  |






Figure 1. Human PLXDC2 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

