

PRODUCT INFORMATION

Target	SPARC
Synonyms	ON; ONT; OI17; BM-40
Description	Recombinant human SPARC Protein with C-terminal 10×His tag
Delivery	In Stock
Uniprot ID	P09486
Expression Host	HEK293
Tag	C-10×His tag
Molecular Characterization	SPARC(Ala18-Ile303) 10×His tag
Molecular Weight	The protein has a predicted molecular mass of 34.1 kDa after removal of the signal peptide. The apparent molecular mass of SPARC-His is approximately 35-55 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	This gene encodes a cysteine-rich acidic matrix-associated protein. The encoded protein is required for the collagen in bone to become calcified but is also involved in extracellular matrix synthesis and promotion of changes to cell shape. The gene product has been associated with tumor suppression but has also been correlated with metastasis based on changes to cell shape which can promote tumor cell invasion. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2015]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Human SPARC Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

