

PRODUCT INFORMATION

| | |
|---|---|
| Tag | C-Flag Tag |
| Target | TAAR6 |
| Synonyms | TA4, TAR4, TAR6, TRAR4, taR-4, taR-6 |
| Description | Human TAAR6 full length protein-synthetic nanodisc |
| Delivery | 6~8weeks |
| Uniprot ID | Q96R18 |
| Expression Host | HEK293 |
| Protein Families | GPCR,Transmembrane,Druggable Genome, |
| Protein Pathways | N/A |
| Molecular Weight | The human full length TAAR6 protein has a MW of 38.5kDa |
| Formulation & Reconstitution | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | This gene encodes a seven-transmembrane G-protein-coupled receptor that likely functions as a receptor for endogenous trace amines. Mutations in this gene may be associated with schizophrenia.[provided by RefSeq, Feb 2010] |
| Usage | Research use only |
| Conjugate | Unconjugated |

