

PRODUCT INFORMATION

| | |
|---|--|
| Target | TPRA1 |
| Synonyms | GPR175, TMEM227, TPRA40 |
| Description | Human TPRA1 full length protein-synthetic nanodisc |
| Delivery | 6~8weeks |
| Uniprot ID | Q86W33 |
| Expression Host | HEK293 |
| Protein Families | Transmembrane,Druggable Genome, |
| Protein Pathways | N/A |
| Molecular Weight | The human full length TPRA1 protein has a MW of 41.1kDa |
| Formulation & Reconstitution | Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Do not use solvents with pH lower than 6.5 in subsequent experiments. |
| Storage & Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Predicted to enable G protein-coupled receptor activity. Predicted to be involved in G protein-coupled receptor signaling pathway. Predicted to act upstream of or within embryonic cleavage and negative regulation of mitotic cell cycle phase transition. Predicted to be integral component of membrane. Predicted to be active in plasma membrane. [provided by Alliance of Genome Resources, Apr 2022] |
| Usage | Research use only |

