

PRODUCT INFORMATION

Target	CD63
Description	The embedded and sliced tumor sections (4µm section) from M-NSG mice implanted with T47-D cells
Host Cells	T47-D
Applications	IHC
Host mice	M-NSG
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Storage & Shipping	Room temperature
Usage	Research use only


SLI100020

Figure 1. IHC analysis of M-NSG T47-D DiSliceX™ SlideSet section (Cat. No. SLI100020) with anti-TROP2 antibody at a concentration of 10 µg/ml. The inset image shows the absence of staining in secondary antibody-only control.



DIMABIO CONFIDENTIAL

