

**PRODUCT INFORMATION**

<b>Target</b>	B7H5, BST1, CCR2, CD33, CD70, CD93, CHI3L1, TREM2
<b>Description</b>	The embedded and sliced tumor sections (4µm section) from M-NSG mice implanted with THP-1 cells
<b>Host Cells</b>	THP-1
<b>Applications</b>	IHC
<b>Host mice</b>	M-NSG
<b>Delivery</b>	1 week
<b>Guarantee</b>	36 months from the date of manufacture
<b>Background</b>	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
<b>Storage &amp; Shipping</b>	Room temperature
<b>Usage</b>	Research use only

DIMABIO CONFIDENTIAL

