

PRODUCT INFORMATION

Target	AFP
Synonyms	Alpha-fetoprotein
Description	Recombinant mouse AFP protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	P02772
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	Mouse AFP(Leu21-Val605) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 66.0 kDa after removal of the signal peptide. The apparent molecular mass of mAFP-His is approximately 55-70 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Predicted to enable fatty acid binding activity and zinc ion binding activity. Acts upstream of or within SMAD protein signal transduction; ovulation from ovarian follicle; and progesterone metabolic process. Located in cytoplasm and extracellular space. Is expressed in several structures, including amniotic fluid; egg cylinder; gut; liver and biliary system; and metanephros. Human ortholog(s) of this gene implicated in hepatitis C. Orthologous to human AFP (alpha fetoprotein). [provided by Alliance of Genome Resources, Apr 2022]
Usage	Research use only
Conjugate	Unconjugated





Figure 1. Mouse AFP Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

