

PRODUCT INFORMATION

Target	CEACAM5
Synonyms	Pregnancy-specific glycoprotein 30;Psg30
Description	Recombinant mouse CEACAM5 protein with C-terminal 6×His tag
Delivery	In Stock
Uniprot ID	Q3UKK2
Expression Host	HEK293
Tag	C-6×His Tag
Molecular Characterization	Mouse CEACAM5(Gln35-Glu947) 6×His tag
Molecular Weight	The protein has a predicted molecular mass of 103.4kDa after removal of the signal peptide. The apparent molecular mass of mCEACAM5-His is approximately 130-180 kDa due to glycosylation.
Purity	The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Cell surface glycoprotein that plays a role in cell adhesion, intracellular signaling and tumor progression. Mediates homophilic and heterophilic cell adhesion with other carcinoembryonic antigen-related cell adhesion molecules, such as CEACAM6. Plays a role as an oncogene by promoting tumor progression; induces resistance to anoikis of colorectal carcinoma cells.[UniProtKB/Swiss-Prot Function]
Usage	Research use only
Conjugate	Unconjugated






Figure 1. Mouse CEACAM5 Protein, His Tag on SDS-PAGE under reducing condition.

DIMABIO CONFIDENTIAL

