

PRODUCT INFORMATION

Uniprot ID	Q13443
Common Name	IMGC 936,IMGC-936,IMGC936
Synonyms	CORD9;MCMP;MDC9;Mltng
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1(K97R,M135Y,S137T,T139E,S325C)
Reactivity	Human
Target	ADAM9
Description	PE-conjugated Anti-ADAM9(izeltabart biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

