

PRODUCT INFORMATION

Clone ID	DM205
Target	ANGPTL3
Synonyms	ANG-5; ANGPT5; ANL3; FHBL2
Host Species	Rabbit
Description	PE-conjugated Anti-ANGPTL3 antibody(DM205); Rabbit mAb
Delivery	Under Development
Uniprot ID	Q9Y5C1
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	This gene encodes a member of a family of secreted proteins that function in angiogenesis. The encoded protein; which is expressed predominantly in the liver; is further processed into an N-terminal coiled-coil domain-containing chain and a C-terminal fibrinogen chain. The N-terminal chain is important for lipid metabolism; while the C-terminal chain may be involved in angiogenesis. Mutations in this gene cause familial hypobetalipoproteinemia type 2.
Usage	Research use only
Conjugate	PE-conjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

