

PRODUCT INFORMATION

Uniprot ID	Q5ZPR3
Common Name	MG-A271, MGA271
Synonyms	B7-H3;CD276;B7 homolog 3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	B7-H3
Description	PE-conjugated Anti-B7-H3 (enoblituzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Research use only□All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.
Usage	
Conjugate	PE-conjugated

