

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	Q7Z7D3
<b>Common Name</b>	AZD-8205; AZD 8205; AZD8205
<b>Synonyms</b>	B7X;B7H4;B7S1;B7-H4;B7h.5;VCTN1;PRO1291
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Host Species</b>	N/A
<b>IgG type</b>	IgG1
<b>Reactivity</b>	Human
<b>Target</b>	B7-H4
<b>Description</b>	PE-conjugated Anti-B7-H4(puxitatug biosimilar) mAb
<b>Delivery</b>	6-8 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

