

PRODUCT INFORMATION

Uniprot ID	Q02223
Common Name	2857916, GSK2857916, J6M0, J6M0-mcMMAF
Synonyms	CD269;TNFRSF17;BCM;BCMA
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	BCMA
Description	PE-conjugated Anti-BCMA (belantamab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

