

PRODUCT INFORMATION

Uniprot ID	Q02223
Common Name	BMS-986352;BMS 986352;BMS986352
Synonyms	BCM;BCMA;CD269;TNFRSF13A□TNFRSF17
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	BCMA
Description	PE-conjugated Anti-BCMA(ispectamab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

