

PRODUCT INFORMATION

Uniprot ID	Q7Z6A9
Common Name	Icatolimab;JS-004;JS004;JS 004;JS-4;JS4;JS 4;TAB-004;TAB004;TAB 004;TAB-4;TAB4;TAB 4
Synonyms	BTLA;BTLA1;CD272
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	BTLA
Description	PE-conjugated Anti-BTLA(tifcemalimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

