

PRODUCT INFORMATION

Uniprot ID	P21730
Common Name	HIB210,HIB-210;HIB 210;MOR044254,MOR-044254;MOR 044254;MOR210;MOR-210;MOR 210;TJ 210;TJ-210;TJ210;TJ210001;WBP2191;WBP-2191;WBP 2191
Synonyms	C5A;C5AR;C5R1;CD88
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	C5AR1
Description	PE-conjugated Anti-C5AR1(izastobart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

