

PRODUCT INFORMATION

Uniprot ID	C7G0W1
Common Name	CAN34D3-65,PF-06443537
Synonyms	IL-31
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Chimeric
IgG type	IgG1
Reactivity	Canis lupus familiaris
Target	IL31
Description	PE-conjugated Anti-Canine IL31(lokivetmab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

