

PRODUCT INFORMATION

Uniprot ID	C7G0W1
Common Name	ZTS-00008183;ZTS00008183;ZTS 00008183;ZTS-8183;ZTS8183;ZTS 8183
Synonyms	Interleukin-31;IL-31
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	N/A
IgG type	IgG2
Reactivity	Canine
Target	Canine IL31
Description	PE-conjugated Anti-Canine IL31(riltovetbart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

