

**PRODUCT INFORMATION**

<b>Uniprot ID</b>	P51681
<b>Common Name</b>	PRO 140, HuPRO-140, PRO-140
<b>Synonyms</b>	CKR5; CCR-5; CD195; CKR-5; CCCKR5; CMKBR5; IDDM22; CC-CKR-5
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	N/A
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG4(P108S)
<b>Reactivity</b>	Human
<b>Target</b>	CCR5
<b>Description</b>	PE-conjugated Anti-CCR5(leronlimab biosimilar) mAb
<b>Delivery</b>	3-4 weeks
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

