

PRODUCT INFORMATION

Uniprot ID	P51685
Common Name	LM108;LM-108;LM 108
Synonyms	CY6;TER1;CCR-8;CKRL1;CDw198;CMKBR8;GPRCY6;CMKBRL2;CC-CKR-8
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CCR8
Description	PE-conjugated Anti-CCR8(cafelkibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

