

PRODUCT INFORMATION

Uniprot ID	P26951
Common Name	IMGN-632;IMGN632;IMGN 632
Synonyms	IL3R;CD123;IL3RX;IL3RY;IL3RAY;hIL-3Ra;IL-3R-alpha
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	N/A
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD123
Description	PE-conjugated Anti-CD123(pivekimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

