

PRODUCT INFORMATION

Uniprot ID	P32927
Common Name	CSL311;CSL-311
Synonyms	CD131;IL3RB;IL5RB;SMDP5;CDw131;betaGMR
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG4
Reactivity	Human
Target	CD131
Description	PE-conjugated Anti-CD131(trabikibart biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

