

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P13726 |
| Common Name | HuMax-TF |
| Synonyms | Tissue factor, TF, F3, Thromboplastin, Coagulation factor III |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | N/A |
| Host Species | Homo sapiens |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD142 |
| Description | PE-conjugated Anti-CD142(tisotumab biosimilar) mAb |
| Delivery | 3-4 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

