

PRODUCT INFORMATION

Uniprot ID	P08571
Common Name	IC14
Synonyms	CD14
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Chimeric
IgG type	IgG4(P1085)
Reactivity	Human
Target	CD14
Description	PE-conjugated Anti-CD14(atibuclimab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

