

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | Q12918 |
| Common Name | IMT009;IMT-009;IMT 009 |
| Synonyms | NKR;CD161;CLEC5B;NKR-P1;NKR-P1A;hNKR-P1A |
| Applications | Flow Cyt |
| Recommended Dilutions | Flow Cyt 1:100 |
| Formulation & Reconstitution | Liquid□PBS with 0.05% Proclin300, 1% BSA |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD161 |
| Description | PE-conjugated Anti-CD161(rezorstobart biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |
| Conjugate | PE-conjugated |

DIMABIO CONFIDENTIAL

