

PRODUCT INFORMATION

Uniprot ID	P20273
Common Name	SM-03;SM03;SM 03;SM-3;SM3;SM 3;SM-06;SM06;SM 06;SM-6;SM6;SM 6
Synonyms	SIGLEC2;SIGLEC-2
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Whole mAb
IgG type	IgG1
Reactivity	Human
Target	CD22
Description	PE-conjugated Anti-CD22(suciraslimab biosimilar) mAb
Delivery	6-8 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

