

PRODUCT INFORMATION

Clone ID	DM63
Target	CD28
Synonyms	CD28; Tp44
Host Species	Rabbit
Description	PE-conjugated Anti-CD28 antibody(DM63); Rabbit mAb
Delivery	3-4 weeks
Uniprot ID	P10747
IgG type	Rabbit IgG
Clonality	Monoclonal
Reactivity	Human
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	The protein encoded by this gene is essential for T-cell proliferation and survival; cytokine production; and T-helper type-2 development. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Usage	Research use only
Conjugate	PE-conjugated

