

PRODUCT INFORMATION

Uniprot ID	P20138
Common Name	CDP-771, CMA-676, WAY-CMA-676,
Synonyms	CD33;SIGLEC3;gp67
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Liquid□PBS with 0.05% Proclin300, 1% BSA
Host Species	Humanized
IgG type	IgG4
Reactivity	Human
Target	CD33
Description	PE-conjugated Anti-CD33 (gemtuzumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

