

PRODUCT INFORMATION

| | |
|---|---|
| Uniprot ID | P20138 |
| Common Name | YMS-5A;YMS5A;YMS 5A |
| Synonyms | p67;SIGLEC3;SIGLEC-3;CD33rSiglec |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | N/A |
| Host Species | Whole mAb |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | CD33 |
| Description | PE-conjugated Anti-CD33(avitotamig biosimilar) mAb |
| Delivery | 6-8 weeks |
| Storage & Shipping | Store at 2°C-8°C for 6 months |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

