

**PRODUCT INFORMATION**

<b>Clone ID</b>	DM135
<b>Target</b>	CD34
<b>Synonyms</b>	CD34
<b>Host Species</b>	Rabbit
<b>Description</b>	PE-conjugated Anti-CD34 antibody(DM135); Rabbit mAb
<b>Delivery</b>	3-4 weeks
<b>Uniprot ID</b>	P28906
<b>IgG type</b>	Rabbit IgG
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	Human
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Purification</b>	Purified from cell culture supernatant by affinity chromatography
<b>Formulation &amp; Reconstitution</b>	Liquid□PBS with 0.05% Proclin300, 1% BSA
<b>Storage &amp; Shipping</b>	Store at 2°C-8°C for 6 months
<b>Background</b>	The protein encoded by this gene may play a role in the attachment of stem cells to the bone marrow extracellular matrix or to stromal cells. This single-pass membrane protein is highly glycosylated and phosphorylated by protein kinase C. Two transcript variants encoding different isoforms have been found for this gene.
<b>Usage</b>	Research use only
<b>Conjugate</b>	PE-conjugated

