

PRODUCT INFORMATION

Uniprot ID	P28907
Common Name	HuMax-CD38
Synonyms	CD38;T10;ADPRC1
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	N/A
Host Species	Homo sapiens
IgG type	IgG1
Reactivity	Human
Target	CD38
Description	PE-conjugated Anti-CD38 (daratumumab biosimilar) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at 2°C-8°C for 6 months
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	PE-conjugated

DIMABIO CONFIDENTIAL

